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MAR

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1. Introduction

■ About MAR

MAR is a fluorescent imaging probe for detection of hypoxia. **MAR** is reduced and makes bright fluorescence under hypoxic condition. **MAR** is able to work under the wide range of pH and to detect mild hypoxia, in which O₂ concentration is around 5%, in living cells.

2. Example of live cell imaging with A549 cells

■ Cell culture

A549 is cultured in DMEM (Dulbecco's modified Eagle's medium) containing 10% fetal bovine serum (FBS) and 1% penicillin streptomycin at 37°C, 5% CO₂ atmosphere.

Seed the cells on to the glass bottom dish coating with poly-L-lysine 2 days before the experiment.

■ Materials Required but not Provided

- Dimethylsulfoxide (DMSO), dehydrated
- PBS

■ Procedure

- ① Dissolve **MAR** in DMSO to 1 mM stock solution.
- ② Dilute the DMSO stock solution with DMEM to 0.5 µM cell stain solution.
- ③ Wash the cell with PBS, add stain solution to the dish and incubate for 1 hour at 37°C, 5% CO₂.
- ④ Wash the cell with PBS, add DMEM to the cell, put the cover slip gently on to the cell and incubate for 3 hours at 37°C, 5% CO₂.
- ⑤ Observe the cells using a fluorescence microscopy.

■ Fluorescent observation

498 nm is suited for excitation wavelength. The wavelength of maximum emission is around 520 nm.

■ Storage

Probes are forwarded under conditions of N₂ atmosphere, dry and frozen state. After receiving, store under -20°C, desiccate and protect from light. We recommend using up DMSO solution of dye.