For Research Use Only **NISPY-3**

Table 1. Product information

Code no.	Product	Contents	Storage	Stability
SK3003 -01	NiSPY-3	1 mg	Freeze-preservation, desiccate and protect from light.	1 year (unopened)

NiSPY-3 (Nitrative Stress Sensing Pyrromethene Dye) is fluorescent reagent for detecting nitrosative stress. **NiSPY-3** reacts with specifically peroxynitrite (ONOO⁻) when it works as detection reagent for reactive oxygen species.

1. About TokyoGreen®-βGlcU(Na)

- NiSPY-3 reacts specifically with peroxynitrite (ONOO⁻) among other reactive oxygen species such as 'OH, O⁻, H₂O₂, ¹O₂, NO etc. Fluorescent intensity is not increased by the existence of hydroxyl radical, singlet oxygen, hydrogen peroxide, hypochlorite, nitric oxide or superoxide.
- Live cell fluorescent imaging is available with NiSPY-3.

2. Principle of the measurement

NiSPY-3 does not have fluorescence in neutral solution. When **NiSPY-3** reacts with peroxynitrite, it becomes to have strong fluorescence (excitation: 490 nm, emission: 515 nm).

3. Contents


NiSPY-3 1mg C₂₃H₁₉BF₂N₄O₄ Mw:464.23

4. Preparation of Reagent

Dissolve all the reagent in 430 μ L of DMSO (5mmol/L), and then dilute 500-5000 times with neutral buffer (final concentration: 1-10 μ M) before use. We recommend to use up at once after probes are diluted.

5. Reference

1. T. Ueno, Y. Urano, H. Kojima and T. Nagano: J. Am. Chem. Soc., 128, 10640-10641 (2006) .



A)Fluorescence spectra of NiSPY-3 solution (10 μ M NiSPY-3 in 0.1 M phosphate buffer pH 7.4 containing 0.1% DMF as a cosolvent) upon addition of peroxynitrite (final 0, 1, 2, 5,10, 20 μ M).

B)Fluorescence response of NiSPY-3 in various ROS generation systems.

"Reprinted with permission from (J. Am. Chem. Soc. 2006, 128, 33, 10640-10641). Copyright (2006) American Chemical Society."

GORYO Chemical, Inc.

EAREE BUILDING 5F 35-100, Nishi-18, Kita-8 Chuo-ku,

Sapporo, 060-0008, JAPAN

URL: http://www.goryochemical.com/