

ProteoGREEN™-gGlu

Table 1. Product information

Catalog no.	Product name	Amount	Storage upon receipt	Stability
GC801	ProteoGREEN-gGlu	20 µg × 10	≤-20°C, keep desiccated and protected from light.	1 year (when unopened and stored as described.)

1. About ProteoGREEN-gGlu

ProteoGREEN-gGlu (gGlu-HMRG) is an activable fluorescent probe to detect expression of γ -glutamyltransferase (GGT). gGlu-HMRG does not emit fluorescence until GGT hydrolyzes it. GGT distributes at cell membrane surfaces and is overexpressed in most cancer cells, so it could be used to detect cancer cells. A green fluorescent product has high membrane permeability and localizes inside the lysosomes




Figure 1. Absorption and emission spectra.
Fluorescence intensity increases 350-fold by GGT

■ Storage

Upon receipt, store the product desiccated and protected from light at $\leq -20^{\circ}\text{C}$. Please use up the solution after dissolving the product. Storing as a solution is not recommended.



■ **Precautions**

Under an acidic environment, ProteoGREEN-gGlu emits fluorescent by the mechanism different from GGT hydrolysis. Therefore, please dilute the probe with a neutral or weak basic solution (pH > 7). The probe also degrades and generates fluorescent products at high temperature, so use the solution promptly after dissolving. To store the solution temporarily, keep the aqueous solution on ice or freeze DMSO solution. Please use the probe rapidly as possible.

2. An example of cell imaging

■ **Preparation of reagent**

1. ProteoGREEN-gGlu is an orange solid. Before opening the vial, spin down the solid to the bottom by a microcentrifuge or by a desktop centrifuge.
2. Warm the vial to the room temperature and add 29.7 μ l of DMSO or PBS to one vial to make 1 mM stock solution. Dilute the stock solution with a neutral buffer or culture medium and use the diluted solution immediately.

■ **Reaction with cells**

1. Dilute 1 mM ProteoGREEN-gGlu with HBSS or RPMI-1640 to 1-2 μ M.
2. Replace the cell culture medium with the cell staining solution and incubate at 37°C, 5% CO₂. for 5-30 minutes. A sufficient increase in fluorescent intensity is observed in GGT-highly expressed cells at 5-15 minutes.
3. Since the background is negligible, it is possible to observe fluorescence without any changes. The fluorescent intensity might increase when enzyme activity proceeds. To prevent such an increase, replace ProteoGREEN-gGlu solution with observation buffer (such as HBSS) and then observe the cells. However, repeated washing would decrease the fluorescent intensity.

■ **Fluorescent observation**

The generated fluorophore can be excited by 488 nm laser and emits fluorescence at 525 nm. For fluorescence microscopy, blue excitation filter cubes for FITC or GFP are suitable.

Table 2. Related Products

Catalog no.	Product name	Major applications
GC6000	HMRG	Hydroxymethyl rhodamine green (green fluorescent probe).
A101-01	MAR	Detection of cellular hypoxia.
GC301	AcidiFluor ORANGE	Detection of lysosomes with orange fluorescence.
SK3001-01	HPF	Detection of hydroxyl radical (\cdot OH) and peroxynitrite (ONOO ⁻).
SK3002-01	APF	Detection of hydroxyl radical (\cdot OH), peroxynitrite (ONOO ⁻) and hypochlorous acid (HOCl).
GC3004-01	OxiORANGE	Detection of hydroxyl radical (\cdot OH) and hypochlorous acid (HOCl).