Behavior of labile ferrous ions and reactive oxygen species during ferroptotic response of cells





Motohiko Sato*, Tasuku Hirayama[¶], Taiga Fujii*, Hideko Nagasawa[¶], Itsushi Minoura* *Goryo Chemical Inc., ¶Laboratory of Pharmaceutical and Medicinal Chemistry, Gifu Pharmaceutical University


info intl@goryochemical.com

Abstract

Ferroptosis is a recently identified iron-dependent regulated cell death (Dixon and Stockwell, 2014). Understanding mechanisms and regulations of this oxidative cell injury may bring us new therapeutic methods for cancer and neurodegenerative diseases. The injury may be caused by reactive oxygen species (ROSs) which are generated by Fenton reaction between labile ferrous ions (Fe2+) and peroxides. Therefore, Fenton reaction should be a key for ferroptotic responses of cells. However, other indirect pathways to generate ROSs should be also considered. To further understand the mechanisms how ROSs are generated during the ferroptotic responses, we tried to visualize labile Fe2+ and ROSs by using cell permeable, activatable fluorescent probes. As previously described, >20 µM of erastin caused ferroptosis of HT-1080 cells. To visualize the responses for an enough time-period (~9 hours), we chose the erastin concentration of 30 µM and observed the behavior of labile Fe2+ and ROSs by fluorescence microscopy. To detect labile Fe2+, both commercially-available chemical probe FeRhoNox-I (RhoNox-I, Hirayama, et al., 2013) and newly-synthesized RhoNox-4 which possesses a higher sensitivity with high specificity to labile Fe2+ were used. To detect ROSs, aminophenyl fluorescein (APF, Setsukinai et al., 2003) which fluoresces by either hydroxyradical, peroxinitrite, or hypochlorous acid and OxiORANGE (Koide et al., 2007) which does by either hydroxyradical or hypochlorous acid were used. The fluorescence signals indicating labile Fe2+ reached the maximum at 3 hours after erasting stimulation. On the other hand, fluorescence signals of APF and OxiORANGE became maximal at 6 hours after the stimulation. The results were consistent with the simple model that hydroxyradical is generated by the reaction between cytoplasmic labile Fe2+ and peroxides and thus causes the ferroptotic cell death.

Dixon and Stockwell 2014, Nat. Chem. Biol. 10:9-17 Hirayama, Okuda, Nagasawa, 2013 Chem. Sci. 4, 1250-1256 Setsukinai, Urano, Kakinuma, Majima, Nagano (2003) J. Biol. Chem. 278:3170-3175 Koide, Urano, Kenmoku, Kojima, Nagano (2007) J. Am. Chem. Soc. 129:10324-10325


Background



precedes the increase in ROSs?


RhoNox-4, a new chemical probe to detect labile Fe²⁺

Specifically reacts with Fe²⁺.




Left, 2 μM of RhoNox-4 was incubated with 1 mM of each ion in 50 mM HEPES buffer. pH=7.4 for 60 min. Right, 2 µM of RhoNox-4 was reacted with 100 µM KO, for 'O, OH generated from 200 µM H₂O₂ and 20 µM FeSO₄, 100 µM NaOCI for HCIO, NO generated by 100 μ M NOC-12, 100 μ M H₂O₃, I mM reduced glutathione, I mM ascorbic acid., and 20 µM FeSO₄. Fluorescence intensities at 575 nm were measured by fluorescence spectrophotometer with 530 nm excitation.

Intracellular Fe²⁺ can be detected Localizes to ER.



HepG2 cells pretreated with/without 100 uM Fe(SO₄)₂(NH₄)₂ for 30 min were washed to remove extracellular Fe2+. Then the cells were reacted with I µM of RhoNox-4 for 30 min. 2,2'-Bipyridine (Bpy, metal-ion chelator) were added with RhoNox-4. The cells were observed by fluorescence microscopy (Ex. 530-560 nm,




HT-1080 cells were co-stained with 125 nM ER-Tracker Green (Molecular Probes, Thermo) and 250 nM RhoNox-4. The cells were observed by fluorescence microscopy (Ex. 460-500 nm, Em. 512-542 nm for ER-Tracker; Ex. 530-560 nm, Em.573-647 nm for RhoNox-4).

Another labile Fe²⁺ probe has been also reported

Allegra T. Aron et al., 2016 J.Am. Chem. Soc., 138: 14338-14346

Chemical probes used



Reactivity of ROS probes

Probes	Alt. name	Ex. (nm)	Em. (nm)	•он	HCIO	ONOO.	H ₂ O ₂
APF	aminophenyl fluorescein	490	515	0	0	0	×
OxiORANGE	MitoHR	553	577	0	0	×	×
HYDROP	NBzF DA	492	518	×	×	×	0

Discussion

- 1) Oxidation of organic molecules such as lipid peroxidation may play more important role than that of H2O2 during the ferroptotic response of HT-1080 cells.
- 2) The decrease in GSH caused by erastin does not directly increase the concentrations of ROSs.

Conclusion

- 1) Chemical probes including RhoNox-4 and ROS probes provide powerful tools to monitor the ferroptotic responses of cells.
- 2) Fe²⁺ concentration increases before the increase in ROSs including H₂O₂.

